

L-Leucin-2-naphthylamid HCl p.A. L-Leucine-2-naphthylamide HCl analytical grade	Kat.Nr./Cat.No. :	27730
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{16}H_{20}N_2O \cdot HCl$
Molekulargewicht molecular weight		292.8
Aussehen appearance		weißes bis hellgraues Pulver white to slightly greyish powder
Reinheit (%) purity	TLC	einheitlich one spot
Löslichkeit solubility	5% in CH_3OH	klar und farblos clear and colorless
Schmelzpunkt (°C) melting point		250 – 260 (Zers./ dec.)
$[\alpha]_D/24^\circ C$	2% in H_2O	+80.0° bis/ to +85.0°
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.