

Malachitgrün Oxalat reinst Malachite green oxalate research grade	Kat.Nr./Cat.No. :	28335
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$(C_{23}H_{25}N_2^+)_2 \cdot C_2O_4^{2-} \cdot 2C_2H_2O_4$
Molekulargewicht molecular weight		927.0
Aussehen appearance		dunkelgrünes Pulver dark green powder
Wasser (%) water	KF	max. 7.0
λ_{max} (nm)	0.0001 % in H ₂ O	618 ± 2
A 1cm/ λ_{max}	0.0001 % in H ₂ O	min. 0.1
TLC	Kieselgel/silicagel 60 F ₂₅₄ n-Butanol/n-butanol 5 Ethanol/ethanol 2 Eisessig/glacial acetic acid 1 Wasser/water 2	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--