

Maltotriose rein pure	Kat.Nr./Cat.No. :	28395
---------------------------------	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{18}H_{32}O_{16}$
Molekulargewicht molecular weight		504.44
Aussehen appearance		weißes krist. Pulver white crystalline powder
Löslichkeit solubility	50 mg/ ml H ₂ O	klare farblose Lösung clear colorless solution
Gehalt (%) assay	HPLC	min. 95.0
[α]²⁵_D	1 % in H ₂ O	+162° bis/ to +166°
Wasser (%) water	KF	max. 6.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.