

α-D-Melezitose reinst	Kat.Nr./Cat.No. :	28550
α-D-Melezitose·research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{18}H_{32}O_{16} \cdot H_2O$
Molekulargewicht molecular weight		522.5
Aussehen appearance		weißes Pulver white powder
Reinheit (%) purity	HPLC	min. 99.0
[α] D/ 20 °C	4 % in H ₂ O	+86.0° bis/ to +90.0°
Wasser (%) water	KF	≤ 5
Schmelzpunkt (°C) melting point		155 – 160
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.