

α-D-Melibiose H₂O reinst	Kat.Nr./Cat.No. :	28560
α-D-Melibiose H ₂ O research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₂ H ₂₂ O ₁₁ ·H ₂ O
Molekulargewicht molecular weight		360.3
Aussehen appearance		weiße Kristalle white crystals
Reinheit purity	TLC	einheitlich one spot
[α]_D/20°C	4% in H ₂ O, 15h	+ 133° bis/ to +137°
Löslichkeit solubility	4% in H ₂ O	klar und farblos clear and colorless
Wasser (%) water	KF	4.5 – 5.5
Schmelzpunkt (°C) melting point		172 – 185
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.