
N,N'-Methylenbisacrylamid 2X p.A. **Kat.Nr./Cat.No. :** **29195**
Bis(acrylamido)methane 2X analytical grade

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₇ H ₁₀ N ₂ O ₂
Molekulargewicht molecular weight		154.2
Aussehen appearance		weiße Nadeln white needles
Gehalt (%) assay	Titration	min. 99.0
A 1cm/ 290 nm	1% in H ₂ O	max. 0.25
pH	1% in H ₂ O	5.5 – 7.5
Löslichkeit solubility	5% in H ₂ O/ 50°C	klar, farblos clear, colorless
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.