

N,N'-Methylenbisacrylamid 2 x Lösung 2% (w/v) N,N'-Methylene bisacrylamide 2 x Lösung 2% (w/v)	Kat.Nr./Cat.No. : 29197
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Beschreibung description	Lösung enthält 2% (w/v) hochreines „Bis“ in deionisiertem Wasser Solution contains 2% (w/v) highly purified „Bis“ in deionized water	
Aussehen appearance		klare farblose Lösung clear colorless solution
A 1cm/ 290 nm	1:1 in H ₂ O	max. 0.25
pH	1:1 in H ₂ O	6.0 – 8.0
Polymerisationszeit (min) polymerizing time	15% Gel, 25°C	max. 30
Lagerung (°C) storage		+4 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.