
Mitomycin C reinst, lyophil. Mitomycin C research grade, lyophil.	Kat.Nr./Cat.No. :	29805
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₅ H ₁₈ N ₄ O ₅
Molekulargewicht molecular weight		334.3
Aussehen appearance		blau-violettes Pulver purple to blue powder
Inhalt content	Jedes Röhrchen enthält 2 mg Mytomycin C und 48 mg NaCl als Verdünnungsmittel. Each vial contains 2 mg Mytomycin C and 48 mg NaCl as diluent.	
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.