

Morpholinoethansulfonsäure • H₂O p.A. Morpholinoethane sulfonic acid monohydrate (MES-monohydrat) analytical grade	Kat.Nr./Cat.No. : 29830
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₆ H ₁₃ NO ₄ S • H ₂ O
Molekulargewicht molecular weight		213.25
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration, calc. on dried substance	min. 99.0
Trocknungsverlust (%) loss on drying		max. 10.0
pH	1% in H ₂ O, 25°C	2.5 – 4.0
A 1cm/260nm	0.1M in H ₂ O	max. 0.05
A 1cm/280nm	0.1M in H ₂ O	max. 0.02
Identität identity	IR	entspricht corresponds
Schwermetalle (ppm) heavy metals	als/ asPb	max. 10
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.