

Morpholinoethansulfonsäure p.A. Morpholinoethane sulfonic acid (MES) analytical grade	Kat.Nr./Cat.No. : 29834
---	--------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₆ H ₁₃ NO ₄ S
Molekulargewicht molecular weight		195.24
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration, getr. Subst. dried subst.	min. 99.0
Trocknungsverlust (%) loss on drying		max. 1.0
pH	10% in H ₂ O	2.5 – 4.0
A 1cm/260nm	10% in H ₂ O	max. 0.1
A 1cm/280nm	10% in H ₂ O	max. 0.1
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10
Identität identity	IR	entspricht corresponds
Ursprung origin		synthetisch synthetic
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.
