

Myoglobin equine lyophil.
research grade

Cat.No.:

29895

Parameter	Method	Specification
Molecular weight		ca. 17 800
Appearance		red lyophilisate
Assay (%)	SDS electrophoresis	min. 95
Solubility	1% in H ₂ O	clear red-brown solution
Water (%)	KF	max. 10
Iron (%)		0.25 – 0.32
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.