

1-Naphthylphosphat-Na-salz p.A.	Kat.Nr./Cat.No. :	30130
1-Naphthyl phosphate-Na-salt analytical grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{10}H_8O_4P \cdot Na \cdot H_2O$
Molekulargewicht molecular weight		264.15
Aussehen appearance		feines weißes Pulver fine white powder
Gehalt (%) assay	Titration	min. 98.0
pH	5% in H_2O	3.5 - 4.5
Wasser (%) water	KF	max. 10.0
Freies Naphthol (%) free Naphthol		max. 0.01
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.