

Natriumcarbonat p.A. Sodium carbonate analytical grade	Kat.Nr./Cat.No. :	30181
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		Na ₂ CO ₃
Molekulargewicht molecular weight		106.0
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	99.5 – 100.5
Identität identity	IR	entspricht corresponds
Chlorid (ppm) chloride	Cl	max. 125
Sulfat (ppm) sulfate	SO ₄	max. 250
Arsen (ppm) arsene	As	max. 5
Eisen (ppm) iron	Fe	max. 50
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 50
Trocknungsverlust (%) loss on drying	300°C	max. 1.0
Lagerung (°C) storage		+15 bis/ to +30

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der Ph.Eur.

The physical/ chemical data corresponds to the requirements of Ph.Eur.

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--