

Natriummolybdat·2H₂O p. A. Sodium molybdate·2H ₂ O analytical grade	Kat.Nr./Cat.No. :	30207
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		Na ₂ MoO ₄ ·2H ₂ O
Molekulargewicht molecular weight		241.95
Aussehen appearance		weißes Pulver oder farblose Kristalle white powder or colorless crystals
Gehalt (%) assay		min. 99.0
Chlorid (ppm)/ chloride	Cl ⁻	max. 50
Eisen (ppm)/ iron	Fe	max. 10
Schwermetalle (ppm) heavy metals	als/ as Pb	max.100
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.