

<b>Natriumthiosulfat-5H<sub>2</sub>O für Silberfärbung</b> p.A. Sodium Thiosulfate Pentahydrate for Silver Staining analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>30218</b>
--	--------------------------	--------------

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> • 5H <sub>2</sub> O
<b>Molekulargewicht</b> molecular weight		248.2
<b>Aussehen</b> appearance		weiße transparente Kristalle white transparent crystals
<b>Gehalt (%)</b> assay	Titration	min. 99.5
<b>pH</b>	10% in H <sub>2</sub> O	6.0 – 8.4
<b>Schwermetalle (ppm)</b> heavy metals	als/ as Pb	max. 5
<b>Kupfer (ppm)</b> Copper	Cu	max. 5
<b>Eisen (ppm)</b> iron	Fe	max. 5
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.