

β-Nicotinamid adenindinucleotidphosphat-Na₂-salz reinst NADP research grade	Kat.Nr./Cat.No. :	30315
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₁ H ₂₆ N ₇ O ₁₇ P ₃ Na ₂
Molekulargewicht molecular weight		787.4
Aussehen appearance		weißes bis leicht gelbliches Pulver white to slightly yellowish powder
NADP·Na₂ (%)	enzym.; 340nm based on dry subst	min. 97.0
NADP (%)	enzym.; 340nm	min. 85.0
NADP (%)	aus/from ε260nm	min. 85.0
Wasser (%) water	K.F.	max. 6.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.
