

β -Nicotinamid adenindinucleotidphosphat-Na₄-salz red. p.A. Kat.Nr./Cat.No. : **30316**
NADPH, analytical grade

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₁ H ₂₆ N ₇ O ₁₇ P ₃ Na ₄
Molekulargewicht molecular weight		833.4
Aussehen appearance		weißes bis leicht hellgelbes Pulver white slightly yellowish powder
NADPH·Na₄ (%)	enzym.; 340nm based on dry subst.	min. 97.0
Assay spectrophotom. (%)	abs. at 260nm ϵ =15 000 based on dry subst.	min. 96.0
NADPH·Na₄ (%)	HPLC	min. 95.0
Wasser (%) water	KF	max. 6.0
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.