

4-Nitrophenyl-β-D-glucuronid reinst 4-Nitrophenyl-β-D-glucuronide research grade	Kat.Nr./Cat.No. : 30725
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₂ H ₁₃ NO ₉
Molekulargewicht molecular weight		315.2
Aussehen appearance		fast weißes Pulver off-white powder
Gehalt (%) assay	HPLC	min. 99.0
[α]²⁰_D	1 % in H ₂ O	−98.0° bis/to −93.0°
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.