

---

<b>4-Nitrophenylphosphat-2CHA-salz p. A.</b>	<b>Kat.Nr./Cat.No. :</b>	<b>30768</b>
4-Nitrophenyl phosphate-2CHA-salt analytical grade		

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		$C_6H_4NO_6P \cdot (C_6H_{14}N)_2$
<b>Molekulargewicht</b> molecular weight		417.44
<b>Aussehen</b> appearance		weiße Kristalle white crystals
<b>Gehalt (%)</b> purity	enzymatisch enzymatic assay	min. 97.0
<b>Freies p-Nitrophenol (%)</b> free p-Nitrophenol		max. 1.0
<b>Identität</b> identity	IR	entspricht corresponds
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.