
Nonactin reinst Nonactin From Streptomyces Griseus, research grade	Kat.Nr./Cat.No. :	30810
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄₀ H ₆₄ O ₁₂
Molekulargewicht molecular weight	für Nonactin for nonactin	736.93
Aussehen appearance		farblose nadelartige Kristalle colourless needle-like crystals
Reinheit (%) purity	TLC Nonactin + Monactin	min. 99
Gehalt Nonactin (%) assay nonactin	TLC	min. 80
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.