

Octyl-β-D-glucopyranosid reinst Octyl-β-D-glucopyranoside research grade	Kat.Nr./Cat.No. :	31055
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₄ H ₂₈ O ₆
Molekulargewicht molecular weight		292.44
Aussehen appearance		weißes kristallines Pulver white crystalline powder
Gehalt (%) assay	HPLC	min. 98.0
[α] D/20°C	5% in H ₂ O	-28° bis/ to -32°
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.