

---

<b>Over Sulfated Chondroitin Sulfate Standard (OSCS)</b>	<b>Kat.Nr./Cat.No. :</b>	<b>31254</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
-------------------------------	--------------------------	---------------------------------------

<b>Identity*</b>	Capillary electrophoresis (CE)	must comply
------------------	--------------------------------	-------------

**For customer information:**

Origin	contaminated heparin
--------	----------------------

Purity	≥ 95 %, based on CE
--------	---------------------

<b>Storage conditions</b>	-15 to -25°C
---------------------------	--------------

**\*Identity assay:**

In CE electropherogram one distinct peak at retention time of 5.2 min must appear.

Capillary: fused silica, 50 µm ID, 67/50 cm

Method: 36 mM monobasic sodium phosphate monohydrate, pH 3.5, wave length 200 nm

Separation conditions: 15 min at 30 KV

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--