

Parthenolide reinst from Chrysanthemum parthenium research grade	Kat.Nr./Cat.No. :	31635
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₅ H ₂₀ O ₃
Molekulargewicht molecular weight		248.3
Aussehen appearance		weiße Kristalle white crystals
Reinheit (%) purity	TLC	min. 98
Löslichkeit solubility	50 mg/ml in CH ₂ Cl ₂	klare, farblose Lösung clear, colorless solution
Schmelzpunkt (°C) melting point		116 – 118
Lagerung (°C) storage		-15 bis/ to -25

Das Produkt ist auch in DMSO (100 mg/ml) und Ethanol (20 mg/ml) löslich.
The product is also soluble in DMSO (100 mg/ml) and Ethanol (20 mg/ml).

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--