

Pentansulfonsäure-Na-salz reinst Pentanesulfonic acid-Na-salt research grade	Kat.Nr./Cat.No. :	31812
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_5H_{11}SO_3 \cdot Na$
Molekulargewicht molecular weight		174.2
Aussehen appearance		weiße Kristallblättchen white crystal leaves
Gehalt (%) assay	Titration	min. 99.0
Identität identity	IR	entspricht corresponds
Löslichkeit solubility	10% in H ₂ O 0.4% in CH ₃ OH	klar und farblos clear and colorless klar und farblos clear and colorless
Wasser (%) water	KF	max. 12.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.