

Phosphoramidon reinst Phosphoramidon research grade	Kat.Nr./Cat.No. :	32753
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{23}H_{32}N_3O_{10}P \cdot Na_2$
Molekulargewicht molecular weight		587.5
Aussehen appearance		weißes bis leicht gelbliches Pulver white to light yellowish powder
Gehalt (%) assay	HPLC	min. 90.0
[α]D/20°C	1% in H ₂ O	- 35° bis/ to - 27°
Trocknungsverlust (%) loss on drying		max. 5.0
Inhibitor Aktivität (µg/ml) inhibitory activity	IC ₅₀	< 0.8
Dünnschichtchromatographie thin-layer chromatography		einheitlich single spot
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.