

o-Phthalaldehyd p.A. o-Phthalaldehyde analytical grade	Kat.Nr./Cat.No. :	32800
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₈ H ₆ O ₂
Molekulargewicht molecular weight		134.14
Aussehen appearance		gelbe Kristalle yellow crystals
Gehalt (%) assay	GC	min. 99.0
Schmelzpunkt (°C) melting point		54 – 57
A 1cm/222nm	0.001% in CH ₃ CN	2.0 ± 0.3
Löslichkeit solubility	15% in C ₂ H ₅ OH	klare gelbe Lösung clear yellow solution
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.