

Progesteron p.A. Progesterone, analytical grade	Kat.Nr./Cat.No. :	33530
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{21}H_{30}O_2$
Molekulargewicht molecular weight		314.45
Aussehen appearance		weißes bis fast weißes feines Pulver white to off white fine powder
Gehalt (%) assay	Ph. Eur.	97.0 – 103.0
$[\alpha]_D/20^\circ C$	Ph. Eur.	+186° bis/ to +194°
verwandte Substanzen related substances	Ph. Eur.	entspricht corresponds
Identität identity	IR	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 0.5
Lagerung (°C) storage		+15 bis/ to +30

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der Ph. Eur.
The physical/ chemical data corresponds to the requirements of Ph. Eur. .

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--