

Proteasome Inhibitor MG-115

Cat.No. :

33767

Parameter	Method	Specification
Amino acid-sequence		Z-Leu-Leu-Nva-CHO
Formula		C ₂₅ H ₃₉ N ₃ O ₅
Molecular weight		461.6
Appearance		white lyophilisate
Purity (%)	HPLC	min. 90
Solubility		soluble in DMSO
Storage/ Stability	powder	at -15 to -25 °C in the dark for more than one year
	solution	at -15 to -25 °C for more than one month if it is protected from moisture and kept in the dark

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.