

Proteasome Substrate

Cat.No. :

33773

Parameter	Method	Specification
Amino acid-sequence		Suc-Leu-Leu-Val-Tyr-MCA
Formula		C ₄₀ H ₅₃ N ₅ O ₁₀
Molecular weight		763.9
Appearance		white lyophilisate
Assay (%)	HPLC	min. 98
Assay	TLC	single spot
Solubility		soluble in DMSO
Storage / Stability	powder	at -15 bis/ to -25°C in the dark for more than one year
	solution	at -15 bis/ to -25°C for more than one month if it is protected from moisture and kept in the dark

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.