

Proteasome Substrate

Cat.No. :

33774

Parameter	Method	Specification
Amino acid-sequence		Z-Leu-Leu-Glu-MCA
Formula		C ₃₅ H ₄₄ N ₄ O ₉
Molecular weight		664.8
Appearance		white powder
Assay (%)	HPLC	min. 98.0
Assay	TLC	single spot
Solubility		soluble in DMSO
Storage / Stability	powder solution	at -15 bis/ to -25°C in the dark for more than one year at -15 bis/ to -25°C for more than one month if it is protected from moisture and kept in the dark

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.