

---

**Pyruvate kinase from rabbit muscle** **Cat.No. : 34085**  
E.C. 2.7.1.40 suspension

---

Parameter	Method	Specification
<b>Molecular weight</b>		ca. 237 000
<b>Appearance</b>		white suspension in 70% sat. ammonium sulfate
<b>Concentration (mg/ml)</b>	Lowry	ca. 20
<b>Activity (U/mg)</b>	Pyruvate kinase	min. 300
<b>Storage (°C)</b>		+2 to +8

**Unit definition**

1 Unit catalyzes the conversion of 1  $\mu$ mol phosphoenolpyruvate to pyruvate per minute at 25°C, pH 7.4.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.