
| | | |
|--|--------------------------|--------------|
| Rapamycin aus Streptomyces hygroscopicus Rapamycin from Streptomyces hygroscopicus, research grade | Kat.Nr./Cat.No. : | 34145 |
|--|--------------------------|--------------|

| Parameter parameter | Methode method | Spezifikation specification |
|---|--------------------------|---|
| Summenformel formula | | C ₅₁ H ₇₉ NO ₁₃ |
| Molekulargewicht molecular weight | | 914.2 |
| Aussehen appearance | | weißes Pulver white powder |
| Gehalt (%) assay | HPLC | min. 97.0 |
| TLC | | min. 97.0 |
| Löslichkeit solubility | 50 mg/ml DMSO | klare farblose Lösung clear colorless solution |
| λ max (nm) | in MeOH | 265, 277, 287 |
| Lagerung (°C) storage | | -15 bis/ to -25 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.