
Rifampicin reinst Rifampicin research grade	Kat.Nr./Cat.No. :	34514
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{43}H_{58}N_4O_{12}$
Molekulargewicht molecular weight		822.96
Aussehen appearance		rotbraunes Pulver red-brownish powder
Gehalt (%) assay		98.0 – 102.0
Identität identity	IR	entspricht corresponds
pH		4.0 – 6.5
Trocknungsverlust (%) loss on drying		max. 1.0
Sulfatasche (%) sulfated ash		max. 0.1
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.