

| | | |
|---|--------------------------|--------------|
| D-Salicin reinst D-Salicin research grade | Kat.Nr./Cat.No. : | 34600 |
|---|--------------------------|--------------|

| Parameter parameter | Methode method | Spezifikation specification |
|---|-------------------------------------|--|
| Summenformel formula | | C ₁₃ H ₁₈ O ₇ |
| Molekulargewicht molecular weight | | 286.3 |
| Aussehen appearance | | weißes feines Pulver white fine powder |
| Gehalt (%) assay | HPLC | min. 99.0 |
| [α]_D/20°C | 5% in H ₂ O | - 59.5° bis/ to -63.5° |
| Schmelzpunkt (°C) melting point | | 200 - 202 |
| Wasser (%) water | Trocknungsverlust loss on drying | max. 0.5 |
| Lagerung (°C) storage | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.