

Serva Blau W Serva Blue W	Kat.Nr./Cat.No. :	35053
-------------------------------------	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{42}H_{46}N_3O_6S_2 \cdot Na$
Molekulargewicht molecular weight		776.0
Aussehen appearance		dunkelviolettes Pulver dark-violet powder
λ_{max} (nm)	0.001 % in Puffer/buffer pH 7	556 - 564
A 1cm/ λ_{max}	0.001% in Puffer/buffer pH 7	min. 0.5
Wasser (%) water	KF	max. 10.0
TLC	Kieselgel/silicagel 60 F254 n-Butanol/n-butanol 5 Ethanol/ethanol 2 Eisessig/glacial acetic acid 1 Wasser/water 2	1 Hauptbande/main spot Rf ca. 0.63 mehrere Nebenbanden/ several secondary spots
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.