

**Serva Violett 17** reinst  
Serva violet 17 research grade

**Kat.Nr./Cat.No.:** 35072

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>41</sub> H <sub>44</sub> N <sub>3</sub> O <sub>6</sub> S <sub>2</sub> ·Na
<b>Molekulargewicht</b> molecular weight		761.9
<b>Aussehen</b> appearance		blauschwarzes Pulver blue-black powder
<b>λmax (nm)</b>	0.001 % in H <sub>2</sub> O	544 – 550
<b>A 1cm/λmax</b>	0.001 % in H <sub>2</sub> O	min. 0.8
<b>Wasser (%)</b> water	KF	max. 10.0
<b>TLC</b>	Kieselgel/silicagel 60 F <sub>254</sub> n-Butanol/ n-butanol Ethanol/ ethanol Eisessig/ glacial acetic acid Wasser/ water	entspricht corresponds
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.