
SERVA DensiStain Blue G Staining Solution, 500 ml 2x Konzentrat	Kat.Nr./Cat.No. :	35078
SERVA DensiStain Blue G Staining Solution, 500 ml 2fold concentrate		

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		blaue kolloidale Lösung blue colloidal solution
$\lambda_{\max.}(\text{nm})$	1:100 in H ₂ O	578 - 586
A 1cm/$\lambda_{\max.}$	1:100 in H ₂ O	min. 0.8
Anwendung application		zum Färben von IEF- und SDS PAGE Gelen for staining of IEF and SDS PAGE gels
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.