
SERVA ProteinStain Fluo-R, Lösung SERVA ProteinStain Fluo-R, solution	Kat.Nr./Cat.No. :	35091
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{72}H_{42}N_6O_{18}S_6 \cdot Na_4Ru$
Molekulargewicht molecular weight		1664.55
Aussehen appearance		rote Lösung red solution
Emissionsmaximum (nm) emission maximum		617
Absorptionsmaximum (nm) absorption maximum		277 437 460
Anwendung application		Proteindetektion in SDS-PAGE protein detection in SDS-PAGE
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.