
Starch (Smithies) hydrolyzed for gel electrophoresis	Kat.Nr./Cat.No. :	35370
--	--------------------------	--------------

Parameter	Method	Specification
Formula		$(C_6H_{10}O_5)_n$
Molecular weight		$(162.14)_n$
Appearance		white powder
Loss on drying (%)	105°C; 3h	max. 20.0
Sulfated ash (%)		max. 1.0
Storage (°C)		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.