
Superoxide dismutase from bovine erythrocytes Min. 3000 units/mg; lyophil., EC 1.15.1.1.	Cat.No. :	35700
----------------------------------------------------------------------------------------------------	------------------	--------------

Parameter	Method	Specification
Molecular weight		ca. 32 000
Appearance		pale blue freeze-dried powder
Solubility	1 mg/ml in 0.05M potassium phosphate buffer, pH 7.8 cont. 0.1mM EDTA	clear and colorless
Activity (U/mg)		min. 3 000
Storage (°C)		-15 to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.