

Taurocholsäure-Na-salz rein Taurocholic acid-Na-salt pure	Kat.Nr./Cat.No. : 35779
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{26}H_{44}NO_7S \cdot Na \cdot n H_2O$
Molekulargewicht molecular weight		537.7 (wasserfrei/ anhydrous)
Aussehen appearance		weißes Pulver white powder
Reinheit (%) purity	TLC	min. 98.0 %
Wasser (%) water	KF	max. 7.0
$[\alpha]_D^{20} / ^\circ C$	3% in H_2O	+ 20.0° bis/ to +24.0°
Cholsäure Natriumsalz (%) cholic acid sodium salt	TLC	max. 1.0
Taurodeoxycholsäure Natriumsalz (%) taurodeoxycholic acid sodium salt	TLC	max. 0.5
Taurin (%) taurine	TLC	max. 1.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.