
Tetrabutylammonium-hydrogensulfat p.A. Tetra-n-butylammonium-hydrogensulfate analytical grade	Kat.Nr./Cat.No. :	35854
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{16}H_{36}N^+ \cdot HSO_4^-$
Molekulargewicht molecular weight		339.5
Aussehen appearance		farblose Kristalle colorless crystals
Gehalt (%) assay	Titration	min. 98.0
Schmelzpunkt (°C) melting point		169 – 172
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.