

<b>N,N,N',N'-Tetramethyl-ethylenediamin (TEMED) reinst</b> N,N,N',N'-Tetramethyl-ethylenediamine research grade	<b>Kat.Nr./Cat.No. :</b>	<b>35925</b>
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>6</sub> H <sub>16</sub> N <sub>2</sub>
<b>Molekulargewicht</b> molecular weight		116.21
<b>Aussehen</b> appearance		klare, farblose bis leicht gelbliche Flüssigkeit clear, colorless to slightly yellow liquid
<b>Gehalt (%)</b> assay	GC	min. 98.5
<b>Wasser (%)</b> water	KF	max. 1.0
<b>nD/20°C</b>		1.4169 - 1.4289
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Bei diesem Produkt kann es gelegentlich zu Trübungen bzw. Ausfällungen kommen, die jedoch die Qualität des Produktes nicht beeinflussen.**

This product may occasionally look cloudy or contain a precipitate which will not influence the quality of the product.

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--