

3,3',5,5'-Tetramethylbenzidin reinst 3,3',5,5'-Tetramethylbenzidine (TMB) research grade	Kat.Nr./Cat.No. :	35926
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₆ H ₂₀ N ₂
Molekulargewicht molecular weight		240.4
Aussehen appearance		hellelfenbeinfarbenes Pulver light ivory-colored powder
Gehalt (%) assay	Titration	98.0 – 102.0
Schmelzpunkt (°C) melting point		168 – 170
Trocknungsverlust (%) loss on drying		max. 0.5
Identität identity	IR	entspricht corresponds
Eisen (ppm) lead	ICP-OES	
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.