
Thrombin from bovine plasma **Cat.No. : 36402**
min. 1000 units/mg protein lyophil. E.C. 3.4.21.5

Parameter	Method	Specification
Molecular weight		ca. 37 000
Appearance		white lyophilisate
Specific Activity (NIH U/mg protein)	Thrombin	min. 1000
Storage (°C)		+2 to +8

Unit definition

One NIH units clots a standard fibrinogen solution in 15 sec. at 37°C.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.