

DL-α-Tocopherol reinst DL- α -Tocopherol research grade	Kat.Nr./Cat.No. : 36570
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₉ H ₅₀ O ₂
Molekulargewicht molecular weight		430.71
Aussehen appearance		klare, farblose bis gelblich braune, viskose, ölige Flüssigkeit clear, colourless or yellowish- brown, viscous, oily liquid
Gehalt (%) assay	GC	96.0 – 102.0
[α]D/20°C	10 % in EtOH, abs.	-0.01° bis/ to +0.01°
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der Ph.Eur.
The physical/ chemical data correspond to the requirements of Ph.Eur.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.