

Trichloressigsäure p.A. Trichloroacetic acid analytical grade	Kat.Nr./Cat.No. :	36910
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂ HCl ₃ O ₂
Molekulargewicht molecular weight		163.40
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration	min. 99.0
Trocknungsverlust (%) loss on drying		max. 0.5
Schmelzpunkt (°C) melting point		min. 57.5
Sulfatasche (%) sulfated ash		max. 0.03
Schwermetalle (ppm) heavy metals		max. 20
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.