
Tris-(2-carboxyethyl)-phosphin-Hydrochlorid, reinst Tris-(2-carboxyethyl)-phosphine-Hydrochlorid- research grade	Kat.Nr./Cat.No. :	36970
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₉ H ₁₆ ClO ₆ P
Molekulargewicht molecular weight		286.70
Aussehen appearance		weiße bis leicht gelbliche Kristalle white to slightly yellowish crystals
Gehalt (%) assay	TLC	min. 99
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.