
Trypsin inhibitor from egg white ca. 20 U/mg
(Ovomucoid), 2 x cryst. lyophil. salt-free

Cat.No. : 37310

Parameter	Method	Specification
Molecular weight		ca. 27 000
Appearance		off white lyophilisate
Solubility	1 mg/ml 0.067M potassium phosphate buffer pH 7.6	clear and colorless
Activity (U/mg)	Trypsin inhibitor	ca. 20
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.