

---

**Trypsin inhibitor from soybean min. 13000 U/mg**  
lyophil.

---

**Cat.No. : 37329**

<b>Parameter</b>	<b>Method</b>	<b>Specification</b>
<b>Molecular weight</b>		ca. 22 000
<b>Appearance</b>		pale yellow lyophilisate
<b>Solubility</b>	1 mg/ml in 0.067 M potassium phosphate buffer pH 7.6	clear and colorless
<b>Activity (U/mg)</b>		min. 13 000
<b>Storage (°C)</b>		+2 to +8

**Unit definition**

1 IU (inhibitor units) inhibits 1 U trypsin as defined by cleavage of 1  $\mu$ mole BAEE (N-benzoyl-L-arginine ethyl ester) per minute.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.